


INSTRUCTIONS

Simple Switch Interface



Pack Contents

Carefully check the contents of the box, which are:

Simple Switch Interface

This manual

Product Description

Simple Switch Interface is a simple, convenient way to connect up to two switches to a computer to access most assistive technology software packages. The unit may be configured to give left click/right click, space/enter, F7/F8 (for Clicker software*) or Gamepad1/2 (for Grid software*).

- USB connectivity.
- No additional driver software required.
- PC and Mac compatibility.
- Supports Left click/ Right click, Space/Enter, F7/F8 and Gamepad1/2.
- Retains mode setting even when powered down.

Compatibility

Simple Switch Interface is compatible with all PC and Mac platforms fitted with a USB socket. It cannot be connected to a PS/2 socket even with an adaptor plug since **Simple Switch Interface** makes use of features which cannot be accessed via the PS/2 connector.

Installation

Ensure your computer is fully turned on and booted up. Locate a spare USB socket on your computer and plug in **Simple Switch Interface**. The unit will be recognised automatically and the relevant drivers loaded. The computer usually tells you when installation is complete, depending upon the operating system you have installed. **Simple Switch Interface** is now ready for use.

Mode Setting

The following table determines the mode of each socket:

Buzzer	Left Socket	Right Socket
1 Beep	Left Click	Right Click
2 Beeps	F7	F8
3 Beeps	Space	Enter
4 Beeps	Gamepad 1	Gamepad 2

Table 1: Socket Settings

To change the mode setting, simply press the small button on the base of the unit (A in Figure 1). The number of beeps emitted by the buzzer indicates the new setting, according to Table 1. The new setting is stored in the unit's memory and is recalled whenever the unit is powered up.





Figure 1

Maintenance

Your Simple Switch Interface has no user serviceable parts. If repair becomes necessary the unit should be returned to Pretorian Technologies or an authorised distributor.

Warranty

Your Simple Switch Interface is warranted against defects in manufacture and component failure. The unit is designed for domestic and educational applications. Use outside of these areas will invalidate the warranty. Unauthorised repair or modification, mechanical abuse, immersion in any liquid or connection to incompatible equipment will also invalidate the warranty.

Pretorian
TECHNOLOGIES

Unit 37 Corringham Road Industrial Estate
Gainsborough Lincolnshire DN21 1QB United Kingdom

Tel +44 (0) 1427 678990 **Fax** +44 (0) 1427 678992

SimplyWorks® is a registered trademark of Pretorian Technologies Ltd.